

V6 Sp2 Enhancements

Module	Item
Clients	Availability of a multi-disciplinary layout on installation
Medicare	New credential functionality for accessing the claiming.com.au portal.

V6 SP2 Production Bug Fix Release Notes

Module	Item Fixed	Issue #
Appointments	Cancelling appointment throws object reference error in V6.0.43356.2	BPA-242
Appointments	Cancelling a meeting in a series causes a crash	BPA-309
Client Details	Selecting an Owner for a Document causes a crash in v6.0.43382.0	BPA-256
Data Maintenance	Removal of Medicare Auth Group entry screen due to changes in authentication functionality.	BPA-301
Invoices	Update GL History record to display Refund created from a payment and vice versa.	BPA-298
Mail Chimp	Mail Chimp updated to work with the new Mail Chimp API and fixed issue with connection not working.	BPA-224
Medicare	Medicare Claim History Report > Clicking "Check Status All" in Medicare claiming report when using grouping causes the report to crash	BPA-254
Medicare	MBS Fee Schedule will only available from the date it is effective not prior to that.	BPA-258
Medicare	Medicare Re-submissions not zeroing previous payment. Update multi-item claims to allow for multiple resubmissions.	BPA-287
Options	Restore Client Layout from File option fixed to display open dialog.	BPA-289
Options	Link to the System Requirements tool in System Options	
Permissions Management	Large sites locking up when resetting permissions - updated permissions to enable one role at a time to be reset.	BPA-184
Referrers	Using the Search in Referrers and then clicking Add New Referrer does not present a blank referrer form	BPA-245
Reminders	SMS Global - End of Life Message and removal of the Register an account with SMS Global link.	BPA-238 BPA-276
Reminders	Sending an adhoc SMS via Bp SMS does not work if previous SMS reply settings were set to Reply to Email	BPA-277
Secure Messaging	Argus email id's not working for all referrers. Check put in place to ensure that the ID saved into Bp Allied is compatible with Argus. Please refer to the point number 1 in the User Interface Details section in the documentation here for further information.	BPA-215
Scheduler	Updated to enable any Email or SMS's sent to update the Client Details> Documents "Added by User" field to be set to "Automated".	BPA-233
User Administration	Google Calendar > Update Google Calendar ID link in Options > User Administration > Google Calendar to the updated documentation link.	BPA-252

V6 Sp2 Documentation Updates

Item
Links to the available videos on Bp Online our online learning centre have been added to the relevant sections in the documentation.

V6 Sp2 Production Known Issues

Module	Item
Medicare	Multi-item Bulk bill, DVA and Patient Claims