



Bp Premier PIP QI Reporting Tool

Change Log

Best Practice Software

This document is proprietary information of Best Practice Software Pty Ltd
Distribution of this information may constitute a breach of our Code of Conduct and may infringe our intellectual property rights.
Document Responsibility: Executive Team, Best Practice Software © Copyright 2016

Document Name Bp Premier PIP QI Change Log

Product Name Bp Premier

Contents

Purpose	3
Change Control	3
Change Log.....	4

1. Purpose

As part of the Practice Incentives Program - Quality Improvement (PIP QI) project, a utility is available for health clinics to enable them to upload aggregated, deidentified patient information to the Federal Government's portal. This document provides detail on all changes that are made to the Practice Incentives Program - Quality Improvement (PIP QI) reporting from within the Bp Premier Reporting Tool.

2. Change Control

Date	Change by	Reason
06/2023	Jamie Atkinson	Initial draft of the document

3. Change Log

Date	Reference	Description
????	--	Initial Release
08/2020	1124	Change the PIP QI report to output JSON instead of XML
09/2020	1185	<ul style="list-style-type: none"> • Measure "denominator" and "numerator" fields: type changed from float to integer. • Metadata "date" field changed to "extraction_datetime". Format for this field defined explicitly in the Description - 2020-08-31 09:17:43.614045 • Metadata "submission_guid" removed. (Not needed because extraction_datetime can perform the same role.) • Minor wording changes to the metadata descriptions. • Remove Disaggregation by Diabetes type in QIM10 (show all diabetes types under the measure under the same values) • Change Practice ID to be BPSnnnnn where nnnnn is site ID • Lowercase measure titles and field titles (e.g. qim01-diabetes, age_group, etc)
09/2020	1204	Patients with ExtractionOptOut set to 1 are now excluded from PIP QI Measures
09/2020	1210	Fixed PIP QI bug where some measures were not being included due to case sensitivity on PHN's end
09/2020	1520	Added SQL indexes to assist with performance for queries used when generating the PIP QI report
02/2022	1786	General PIP QI Reporting Tool fixes