

Create a Custom Preparation



Create a custom preparation

To create a custom preparation from the main screen of Bp Premier:

1. Select **Setup > Custom preparations**.
The **Custom prescription preparations** screen will appear.
2. To create a custom preparation:
 - i. Click **Add** to create a new preparation; or
 - ii. Select **Item name** and click **Clone** to copy and modify an existing preparation.

The **Custom prescription item** screen will appear.

3. Enter a **Name** and **Strength** for the preparation.
4. Select the **Form** of the preparation (e.g., Scalp, skin treatment).
5. Click **Add** to add an ingredient from the Ingredient screen.
6. Enter a **Strength** and **Unit** for each ingredient and click **Save** to save the ingredient to the Ingredients list.
7. Repeat **Steps 5-6** to add as many ingredients as required.
8. Enter **Other details** that patients and providers need to know about this preparation.
9. Select the **Route** of administration (e.g., Topical).
10. Tick boxes to check for **Drug interactions** and **Disease interactions** when this Rx is added to a patient's record.
11. Enter the **Quantity** and **Units** for a prescription of this preparation, and the default number of **Repeats**.
12. Select the **Restriction** that applies to this preparation (e.g., No PBS).
13. Click **Save**.

Custom prescription item

Name: Scalp lotion

Strength: Mild

Form: Scalp, skin treatment

Ingredients:

Ingredient name	Strength	Unit	Add
Acacia, powdered	1.0	1.0	Edit
Isopropyl lanolate	1.0	1.0	Delete

Other details:

Route: Topical

Check for: Drug interactions Disease interactions

Quantity: 1 Unit: 100g

Repeats: 0

Restriction: Non PBS

Save Cancel

Prescribe a custom preparation

To prescribe a custom preparation from a patient's record:

1. Select **Current Rx** from the tree on the left-hand side and click **Add**.
2. In the New Rx screen, select **Custom** from the radio buttons along the top.
3. Select the **custom preparation** from the list.
4. Continue [completing the new Rx](#) like any other prescription.

New Rx

Product name Therapeutic class Custom Favourites

Search for: [] Available formulations: Show PBS/RPBS listed only

Drug name	Product name	Quantity	Rpts	Restriction	BPP	TGP/SPC
Scalp lotion	Mils Scalp, skin treatment	1x100g	0	Non PBS/OTC	\$ 0.00	\$ 0.00

Scalp lotion Scalp, skin treatment

Generic name:
Acacia, powdered 1.01.0
Isopropyl lanolate 1.01.0

Schedule:
PBS Listing: Non PBS - 1x100g and 0 repeats

Add to favourites

No allergy/adverse drug reaction information has been recorded!

Same drug class Equivalent products Product Information
NPS RADAR Allergies/Reactions CMI

< Back Next > Cancel