

Bp VIP.net Topaz SP1 Rev 1 Release Notes

These release notes describe all enhancements and modifications made to Bp VIP.net since version Topaz SP1 (2.1.529.008).

Release Date	20 June 2022
Release version	Topaz SP1 Revision 1 (version 2.1.529.009)
Notes last updated	20 June 2022
What is in this release?	Bp VIP.net 2.1.529.009 contains Day Stay updates for Australian practices only.
Which version can I upgrade from?	There is no version dependency for upgrading to 2.1.529.009. If you experience difficulties upgrading from your version, contact Best Practice Software Specialist Product Support. The operating system and database requirements have not changed since the last release.
When should I upgrade?	Best Practice Software recommend running the program upgrade outside of business hours to minimise the impact of the upgrade on your practice.
How do I upgrade?	Upgrade instructions are included in these release notes.
What do I do after upgrading?	You do not need to perform any configuration after you upgrade.

Copyright Statement

This document is classified as commercial-in-confidence. Unauthorised distribution of this information may constitute a breach of our Code of Conduct, and may infringe our intellectual property rights. This information is collected and managed in accordance with our Privacy Policy, available on our [website](#).

© Copyright 2022

Upgrade Bp VIP.net to Topaz Revision 1 (2.1.529.009)

You should be able to upgrade to 2.1.529.009 from any version of Bp VIP.net. Best Practice Software recommend that the upgrade be performed outside of normal business hours.

Notes on upgrading

- If you did not upgrade to the most recent version Topaz SP1 2.1.529.008, [review the release notes](#) to see the changes in that version.
- Copy the install file to and run the upgrade from the **SQL Server database computer** first. This will upgrade the Bp VIP.net database structure in preparation for terminal server or workstation upgrades.
- If your practice uses Terminal Services, you must run the installation from the Terminal Server (after upgrading the database server).
- Log in to the Bp VIP.net server computer as a Windows **administrator**. You may require Windows administrator permissions to install and update some Bp VIP.net components.
- If you are upgrading from a version prior to **525.018**, ensure you have updated your MIMS data to the latest release, to enable Active Ingredient Prescribing activation. If you are already using Active Ingredient Prescribing, no action is necessary.

Upgrade to version 2.1.529.009

1. Check the version of Microsoft .NET Framework on the Bp VIP.net server and workstations. Topaz SP1 Rev 1 requires a minimum .NET Framework version of 4.8.
2. **Back up your existing database** before upgrading and **test that the backup restores successfully**. Information on how to back up and restore your database is available on the knowledge base.
3. Download [the setup.zip file here](#). Your browser may warn that the file is a security risk. You can safely keep this file.
4. **Rename** the file from Setup.2.1.529.009.zip to **setup.zip**.
5. Copy the setup.zip file to the folder nominated in **Setup > Facility Preferences > General Preferences > Login/Tasks** tab > **New version location** field. This field must contain a UNC path to a network location accessible by Bp VIP.net workstations.
6. Log in to the database server as an administrator user. Bp VIP.net will detect that a new setup.zip file exists and begin upgrading the server and database.
7. When a user next logs in to a workstation, Bp VIP.net will detect that a new setup.zip file exists in this folder and will automatically upgrade the version of Bp VIP.net on the workstation.
8. If your practice uses a MIMS database, check your version of MIMS and update to the latest version if available.

Day Stay Reporting

Day stay reporting requirements have been updated in compliance with Australian state standards.

ICD10 Version 12 Download

The ICD10 Version 12 import file is supplied by Best Practice Software support. It will be given in a .zip file format, and must be unzipped once saved to the desktop. The file within the unzipped folder to import into Bp VIP.net. is **package.csv**.

To import the file:

- Go to **File > Utilities > Live Update**.
- Select **Update > Coding > Browse**.
- Select the ellipses button [...] to browse for the previously saved file, **package.csv** and **Open**.
- Select **Next**.
- Tick the ICD10 Version 12 files to update.
- Select **Next**.

Note: *The ICD10 Version 12 update will take 5-6 minutes to process and should not be interrupted.*
